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
1. Air block applied during the final polymerisation

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
3. Macro-geography can be corrected by using fine grain diamond burs (red ring: 861514012; 853514008)...

4. ... or a multi-blade bur (CS9904).

5. Diamond and paper strips and mini-discs for interproximal area

6. Diamond silicon rubber (Shiny 14) for pre-polishing

7. Micro-geography can be checked by marking with pencil

8. ... and is enhanced with a carbundum bur (Shiny 21) or a multi-blade (CS9904)

9. A 3 micron diamond paste (Shiny A) is applied ...

10. ... with a soft hair goat brush (Shiny S)

11. A 1 micron diamond paste (Shiny B) is applied ...

12. ... with a soft hair goat brush (Shiny S)

13. Hi-polishing with Al-Oxide (Shiny C)

14. ...and a felt (Shiny F)

15. Only convex areas (marked with stars) need to be polished with aluminium oxide paste (Shiny C)

16. Final polished composite restoration

17. Impregnated brushes (silicon carbide)

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(GB) FINISHING KIT, ENAMEL PLUS SHINY

A modern aesthetic composite restoration foresees the accurate preparation of the cavity, the use of a stratification technique through dentin fluorescent masses and opalescent enamel masses, and an accurate finishing and polishing phase, able to give the restoration surface a similar feature to that of a natural tooth.

Rotative instruments

The rotative instruments used in the finishing and polishing of composite restorations consist of diamond burs and rubbers. Diamond burs have different grain measures and are usually used in sequences: those with a bigger grain are used for the carving and first working of the restoration, those with a smaller grain, for the finishing. Rubbers can be used exclusively in the last phase of finishing, being aware not to push too much on the tooth surface, in order not to overheat it, creating pulpar damage, and in order not to burn the resin matrix of the composite material.

Abrasive discs and strips

The clinic use of flexible discs is limited because they are not able to reach certain anatomic areas, such as ditches, and they generally reduce characterisations, during the carving phase of restoration, of hardly achieved surfaces. Their use must therefore be limited to incisal and interproximal areas. Abrasive strips are used for the finishing and polishing of interproximal edges. Great care must be taken during their use, in order not to damage the parodontal legament and the gingival papilla, and in order to respect the holding of an interproximal contact point. These are used in sequences as well, and foresee different passages to obtain a good final result.

Polishing pastes

Hybrid composites of the last generation respond to polishing very well, through diamond pastes and aluminium oxide pastes. The Shiny pastes are available in two different grain measures: 3 microns (paste A), and 1 micron (paste B), to be used in sequence with brushes. A final polishing is achieved with felts, using the aluminium oxide paste (paste C).

CLINICAL PROCEDURE

Fig. 1-2. Before starting the finishing of a composite reconstruction, it is necessary to apply an Air Block (Shiny G), when restoration is finished and before the final light curing. This glycerin-based product eliminates the oxygen from the surface (that inhibits the complete curing of the composite), and increases the resistance of the surface material. Tests have shown that using an air block in combination with diamond pastes, the surface hardness is increasing of 10% and its deviation from medium value decreases from 1,2% to 0,8%, thus the surface is more even. (Fig. 2 dev=hardness deviation from medium value, AB= Air Block, DP= polishing with diamond paste).

Fig. 3-4. The anatomy of an anterior tooth is often characterised by a rich surface morphology, represented by ditches and lumps (macrogeography), and by growth lines of the enamel (microgeography). The first finishing phase of a composite restoration foresees the use of Ena Shiny diamond burs in a fine grain and/or multiblade burs. With these burs, used at a low speed, a revision of the surface morphology, obtained through carving, is made; the

final macrogeography is a result. In order to facilitate this phase, the obtained surface can be evidenced, during carving and finishing (macrogeography), by using the Shiny L varnish. The colour subtraction technique allows to clearly visualize the depth of ditches and the disposition of lumps.

Fig. 5. In interproximal areas, a started hand job can be beautifully finished reducing, in areas difficult to access, the build up of bacterial plaque, thanks to diamond strips (available in fine grain red and extra-fine yellow, in 2 and 4 mm height) or paper abrasive strips (Shiny SM and SMS, with double grains, medium-coarse, with two widths, Shiny SF and SFS with double grains fine and super-fine, with two widths) and discs (Shiny DM medium and DF super-fine).

Fig. 6. Finishing is ultimated by using a diamond silicon rubber that has the shape of a flame (Shiny 14), dunked in pumice.

This instrument, used at a low speed and under a spray of water, levigates the surface without erasing the macrogeography.

Fig. 7-8 The microgeography is created with the help of a Shiny carborundum bur (Shiny 21) o multi-blade bur (CS9904), which should be used with great care. By rubbing the dental surface with an articulation paper, the intensity of the created macrogeography is evidenced.

Fig. 9-10-11-12. The pre-finishing, obtained with the silicon rubber, prepares the restoration to a proper finishing, that renders the surface shiny, without destroying the micro and macrogeography. For this, a goat brush (Shiny S) is used, together with diamond pastes Shiny A and B, with different grains (3 microns and 1 micron). The surface must first be treated without water, and at a very low speed, with short working times, following the inclination illustrated in the picture, in order to avoid over-heat; afterwards, a heavy spray of water may be used, increasing speed slightly. For posterior, Shiny 41-42-43 brushes impregnated with silicone carbide can be used

Fig. 13-14-15. Convex areas of the surface, more exposed to using, always appear shinier than the concave ones; therefore, with a round felt (Shiny F), or a disc felt (Shiny FD), and a paste Shiny C in aluminium oxide, we obtain a brighter reflected light in these areas. The paste can be used every time a patient returns for prophylaxis, to polish composite restorations (maintaining phase).

Fig. 16. These different working phases that can be done in one sitting, if made properly, allow the clinician to obtain a restoration that really simulates the real tooth, with a guarantee of the respect of mechanical characteristics of used material, of dental tissues and parodontal tissues.

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Precautions:

- Use rotating instruments at low speed and intermittently. Use in conjunction with water spray to reduce heat production.
- Automatic: in ultrasonic cleaner min. 15 minutes. Manual: cleaning agent validated for burs, disposable tissue, running water.
- Diamond strips, brushes and rubbers can are reusable. Sterilize before each use in auto-

clave at 134°C/273°F min. 3 minutes. The other instruments are disposable.

- Use spray with suction to remove glycerin and varnish.

Side effects and contra-indications

- Do not use discs on the vestibular surfaces in order to avoid flattening the surface (they can be used interproximally and on the incisal edges).
- Use rotating instruments at low speed to avoid overheating the pulp (or the tooth structure - risk of necrosis) and the composite (risk to burn the resin matrix).

Take care:

- SHINY L: Shake well before use. Close cap tightly after use.
- A-B-C-G PASTE: If the material has not been used for a long time, dispose of the dried-out material at the tip of the syringe. Close cap tightly after use. MSDS available on web site www.micerium.com

Note: instruments and materials do not contain hazardous components. They are safe to use when instructions are followed.

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File: Shiny Istr 7 lingue EU V2_2-2012 © Micerium S.p.A